

PALPITATIONS

ARRYTHMIAS

CHEST PAINS

THE CHILTERN BOX

Portable ECG Recorder



The device consists of
a tape recorder box with two leads
attached to a jack-plug

At the onset of symptoms
the patient attaches the two
electrode clips to either an arm and a leg
or both arms—

and then the patient pushes the
jack plug into the Box which
automatically starts the 15 minute
recording of the ECG

Analysing the Recording

The tape recorder in the Chiltern Box is used to play the tape back through a de-coder unit connected to any ECG monitor.

The decoder unit has 3 standard ECG type contacts for connection to the monitor—these are duplicated to permit simultaneous connection to a standard ECG machine so that permanent tracings can be obtained of the 15 minute recording.

Consider the Advantages—


SAVES PHYSICIAN TIME
SAVES TECHNICIAN TIME
LOW INITIAL COST
LOW RUNNING COST

Only one test is performed—and analysed
Compatible with all standard ECG equipment
Standard tapes and batteries have long lives

Ask for full details from:

CHILTERN INTERNATIONAL LIMITED

'Kynance', Manor Road, Penn, Buckinghamshire, HP10 8JB, England
Telephone: Penn (049481) 2177




... a different approach to Patient Monitoring

Electronic monitoring is now an accepted part of patient care. ECG is the most regularly monitored signal and the table below approximates to the percentage of patients monitored electronically for the more common parameters.

Today the cost of quality medical instrumentation does not allow there to be equipment in the intensive or

Coronary Care Unit which is not fully used or is unavailable for use on one patient because it is locked into a system on another. The S&W 8000 System allows total flexibility to match the instrumentation to the patient's need with the added advantages of: ■ Progressive expansion of the monitoring system as new modules for new parameters are developed ■ Modules may be used with a low cost power supply as monitors in their own right; for example, pressure, temperature, respiration ■ May be used in any department in the hospital; C.C.U., I.T.U., Operating Theatres, Neo-natal units ■ Ease of servicing in event of unsatisfactory performances.



S&W

Simonsen & Weel Limited, Unit 1, Ruxley Corner, Sidcup, Kent DA14 5BL. Tel: 01-309 0433.